

SKY pro guide



SIMPLE

SAFE


FAST

Guided and highly precise with just one sleeve.

SKY pro guide

2

Your introduction to precision in guided implantology



SKY pro guide surgical tray (instruments included)

REF SKYPROGST1


Instruments for all SKY implants
with Ø 3.5, Ø 4.0 and Ø 4.5.



SKY guided fixation pin
Ø 1.35 mm L 27mm
REF GFPL27D135





SKY guided fixation pin sleeve
Ø 1.35 mm L 6mm
REF GFSL06D135



The fixation pins enable secure fixation of the drilling template.

Principle of function

Maximum safety when reaching the drilling depth and angulation by guiding the drills through the high-precision drill sleeve.



Function description




Practical: the lower part of the tray contains ratchets and tools for inserting the implants.



One system for all SKY implants from Ø 3.5 – 4.5 mm. In addition, smaller or larger individual bore diameters, such as for pilot holes, can also be achieved in a guided manner.

Case example: bSKY4514



SKY pro guide

Component overview

SKY
PRO GUIDE

SKY pilot drills, guided



PCDXLXD35n PCDXLXD400 PCDXLXD450

SKY twist drills, guided



TDXL52D225 TDXL08D225 TDXL10D225 TDXL12D225 TDXL14D225 TDXL16D225

SKY final drills, guided – hard bone D1



colour code

3.5 N 3.5
4.0
4.5

L 5.2 mm

FD1L52D350
FD1L52D400
FD1L52D450

L 8 mm

FD1L08D350
FD1L08D400
FD1L08D450

L 10 mm

FD1L10D350
FD1L10D400
FD1L10D450

L 12 mm

FD1L12D350
FD1L12D400
FD1L12D450

L 14 mm

FD1L14D350
FD1L14D400
FD1L14D450

L 16 mm

FD1L16D350
FD1L16D400
—

SKY final drills, guided – medium hard / soft bone D2, D3, D4



colour code

3.5 N 3.5
4.0
4.5

L 5.2 mm

FD2L52D350
FD2L52D400
FD2L52D450

L 8 mm

FD2L08D350
FD2L08D400
FD2L08D450

L 10 mm

FD2L10D350
FD2L10D400
FD2L10D450

L 12 mm

FD2L12D350
FD2L12D400
FD2L12D450


L 14 mm

FD2L14D350
FD2L14D400
FD2L14D450

L 16 mm

FD2L16D350
FD2L16D400
—

SKY crestal drills, guided



CRDXLXD35n CRDXLXD350 CRDXLXD400 CRDXLXD450

Implant mounters, guided – copaSKY



COPMLXD475 COPMLXDSET COPAIMSM163

Implant mounters, guided – SKY



SKYMLXD475 SKYMLXDSET SKY-DS04

Implant mounters



SKY-WTK1 SKY-WTK6 SKY-STK1 COPACTK6

Prosthetic keys




SKYC-WM6 SKY-SD22 SKY-SD28

Drill sleeves – pro guide



GDLSL06D475 GDLSL06D47510

SKY pro guide – surgical tray



instruments
not included
SKYPROGST0

instruments
included
SKYPROGST1

SKY Torque Wrench Pro



SKYTWPRO SKYTWCON

Fixation pin sleeves




GFSL06D135 GFSL06D13510

Twist drill for fixation pins, guided



SKYDT13L

Fixation pins



GFPL27D135 GFPL27D1354

Orders are subject to the General Terms and Conditions of bredent medical GmbH & Co. KG: <https://www.bredent-implants.com/general-terms-of-business-gtb/>
Illustrations not to scale. All prices in € plus VAT, plus shipping if applicable.

0091039GB-20200714 Mistake and subject to change reserved



DENTAL INNOVATIONS
SINCE 1974

bredent
group

bredent medical GmbH & Co. KG · Weissenhorner Str. 2 · 89250 Senden · Germany · T: +49 7309 872-441 · F: +49 7309 872-444 · www.bredent-medical.com · @: info-medical@bredent.com